	Standard Operating Procedure	SOP Number D-1020	Revision 0
	Determination of Apigenin by HPLC-UV	Effective Date 08/16/23	Page Page 1 of 7
Written by/ Date CJF 08-15-23	Reviewed by/ Date SAS 08/15/23	Approved by/ Date SSS 08/16/23	
Title: Analytical Development Scientist	Title: Analytical Development Scientist	Title: Quality Control Director	

## 1.0 Purpose

This document describes the analytical procedure for the determination of Apigenin in raw materials and finished products.

## 2.0 Scope

This procedure applies to the identification and quantification of Apigenin in raw materials and finished products. This method was validated under protocol PRTCL-23-0050.

## 3.0 Responsibility

- 3.1 It is the responsibility of QC and Analytical chemists who have verified their ability to execute this procedure to follow this procedure.
- 3.2 It is the responsibility of QC Laboratory Management to implement this procedure and to ensure that the procedure is being followed.
- 3.3 It is the responsibility of QC Laboratory Management and/or Analytical Development Personnel to keep this procedure current with the associated monographs and laboratory practices.

## 4.0 Definitions

- 4.1 **QC** – Quality Control
- 4.2 **MeOH** – Methanol
- 4.3 **DMSO** – Dimethylsulfoxide

- 4.4 **H<sub>3</sub>PO<sub>4</sub>** – Phosphoric Acid
- 4.5 **ACS** – American Chemical Society
- 4.6 **HPLC** – High Performance Liquid Chromatography
- 4.7 **UV-Vis** – Ultraviolet-Visible (Detection)
- 4.8 **SPP** – Superficially Porous Particle (Solid Core Column Stationary Phase)

## 5.0 References

- 5.1 PRTCL-23-0050, Protocol, Validation of an Analytical Method for the Determination of Apigenin by HPLC-UV

## 6.0 Supplies

- 6.1 Chemicals – All reagents are ACS grade or better
  - 6.1.1 MeOH
  - 6.1.2 DMSO
  - 6.1.3 H<sub>3</sub>PO<sub>4</sub>
  - 6.1.4 Apigenin Reference Standard
- 6.2 Supplies and Glassware
  - 6.2.1 HPLC vials, 12mm X 32mm with screw cap enclosures w/ septa
  - 6.2.2 Volumetric glassware
  - 6.2.3 Weigh paper
  - 6.2.4 Syringes with 0.45μ nylon syringe filters

### 6.3 Equipment

6.3.1 Suitable gradient HPLC system consisting of a pump, autosampler, column compartment and UV-Vis detector with a chromatographic data handling system

6.3.2 Analytical Balance

6.3.3 Micro Analytical Balance

6.3.4 Wrist Action Shaker

## 7.0 Procedure

### 7.1 Mobile Phase, Extraction Solvent & Diluent Preparation

#### 7.1.1 Mobile Phase

7.1.1.1 Mobile Phase A: Add 1000 $\mu$ l H<sub>3</sub>PO<sub>4</sub> to 1000ml Milli-Q water and mix well.

7.1.1.2 Mobile Phase B: Use MeOH.

#### 7.1.2 Extraction Solvent

7.1.2.1 Add equal volumes of MeOH and DMSO and mix well.

#### 7.1.3 Diluent

7.1.3.1 Use MeOH.

### 7.2 Standard Preparation

7.2.1 Prepare standard stock at ~0.05 mg/mL analyte in Extraction Solvent. Mechanically shake stock for 15 minutes at ~ 1/2 volume then QS to volume. Shake vigorously then dilute stock 1:10 with Diluent. Shake vigorously then filter a 5mL aliquot for analysis, discarding the first 3-4mL of filtrate.

7.3 Sample Preparation

7.3.1 Specific sample testing details are provided in each product profile. If a specific testing details section is not available, follow preparation procedure as described below, maintaining concentration within the linear range of this method.

7.3.2 The validated linear range for the analytical method is 0.0030 – 0.0079 mg/mL.

7.3.3 Pool at least 10 dosage units and homogenize as appropriate (e.g. grind tablets / capsule fill / powders / stick pack contents by mortar and pestle, cryogenically powder and dissolve gummies, etc.) Extract sufficient sample (based on the raw material manufacturer assay value / finished product profile) with Extraction Solvent in order to generate a sample stock that is ~0.05 mg/mL analyte. Mechanically shake stock for 15 minutes at ~ ½ volume then QS to volume. Shake vigorously then dilute stock 1:10 with Diluent. Shake vigorously then filter a 5mL aliquot for analysis, discarding the first 3-4mL of filtrate.

7.4 HPLC Parameters

7.4.1 Column: Waters Cortecs Shield RP18, 2.7µm (SPP), 4.6 x 150mm

7.4.2 Column Temperature: 40°C

7.4.3 Flow rate: 0.7 mL/min

7.4.4 Run Time: 15 minutes

7.4.5 Mobile Phase Gradient:

7.4.6 Wavelength: 339 nm

7.4.7 Injection Volume: 5 µL

Time, min	% A	% B
0	55	45
3	55	45
10	15	85
10.1	55	45
15	55	45

7.4.8 Suggested 3-D Spectral Range (for Identification): 210nm - 450nm

7.5 Recommended Sequence

- 7.5.1 Make at least 2 injections of the Diluent.
- 7.5.2 Make at least five (5) injections of Working Standard.
- 7.5.3 Make a single injection of each Sample Preparation.
- 7.5.4 Make a single injection of the Working Standard after every ten (10) sample injections or at the end of the run.

7.6 System Suitability Requirements

- 7.6.1 The %RSD of five (5) consecutive standard injections is NMT 2%.
- 7.6.2 The %RSD of all standard injections is NMT 3%.

7.7 Example calculations for determining % Assay:

$$7.7.1 \quad \% = \frac{R_u}{R_s} \times \frac{Wt_{std} \times P}{V_{std}} \times \frac{V_{spl}}{SA} \times 100 \times SS/LA$$


- $R_u$  Sample peak area
- $R_s$  Mean (All) standard peak area
- $Wt_{std}$  Weight of the reference standard, mg
- $V_{std}$  Volume of the standard preparation accounting for dilutions, mL
- $P$  Purity of the reference standard in decimal format
- $SA$  Sample amount, mg
- $V_{spl}$  Volume of the sample preparation accounting for dilutions, mL
- $SS$  Serving Size, mg (Use "1" for raw materials.)
- $LA$  Label Amount, mg (Use "1" for raw materials.)

7.8 System Wash, Column Wash and Column Storage


7.8.1 Wash system and column with 90% ACN. Store column on 50% ACN.

## 8.0 Chromatograms


### 8.1 Typical Diluent Chromatogram



### 8.2 Typical Standard Chromatogram



### 8.3 Typical Finished Product Chromatogram



### 9.0 Revision History

Revision	Date	Description of Changes	CCR #	By
0	08/08/23	New procedure.	N/A	C. Perry