	Standard Operating Procedure		SOP Number D-1024	Revision 1
	Determination of Mannose by HPLC/UV		Effective Date	Page Page 1 of 11
Written by/ Date <i>MMA 10/09/25</i>		Reviewed by/ Date <i>ATS 10/09/25</i>		Approved by/ Date <i>Franz Roly 21-Oct-2025</i>
Title: QC Microbiologist II		Title: QC Laboratory Manager		Title: QA/QC Director ^①

^① Franz Roly, HBI Amersale QC Manager approving in absence of Nicole Motley. ~~FR~~ 21-Oct-2025

1.0 Purpose

The purpose of this procedure is to define the method for the determination of mannose in raw materials and finished products by HPLC/UV.

2.0 Scope

This procedure applies to the determination of mannose by HPLC/UV in the QC laboratory at HBI Ion Labs.

3.0 Responsibility

- 3.1 It is the responsibility of QC Chemists to follow this procedure.
- 3.2 It is the responsibility of QC Laboratory Management to ensure that this procedure is being followed.
- 3.3 It is the responsibility of QC Laboratory Management and/or Analytical Development to keep this procedure aligned with current practices.

4.0 Definitions

- 4.1 **HPLC/UV** – High Performance Liquid Chromatography with Ultraviolet Detection
- 4.2 **QC** – Quality control
- 4.3 **PMP** – 3-methyl-1-phenyl-2-pyrazoline-5-one
- 4.4 **HCl** – Hydrochloric acid

Standard Operating Procedure Determination of Mannose by HPLC/UV	SOP No D-1024	Rev 1	Page 2 of 11
--	--------------------------------	------------------------	-------------------------------

- 4.5 **NaOH** – Sodium hydroxide
- 4.6 **H₃PO₄** – Phosphoric acid (~85%)
- 4.7 **K₂HPO₄** – Potassium phosphate dibasic
- 4.8 **H₂O** – Deionized water (>18MΩ·cm)
- 4.9 **DAD** – Diode array detector
- 4.10 **ACN** – Acetonitrile
- 4.11 **MeOH** - Methanol

5.0 References

- 5.1 PRTCL-24-0015, Protocol, Validation of a Method for the Determination of Mannose by HPLC/UV.
- 5.2 Gonzalez NM, Fitch A., Al-Bazi J. (2020) Development of a RP-HPLC method for the determination of glucose in *Shewanella oneidensis* cultures utilizing 1-phenyl-3-methyl-5-pyrazolone derivatization, PLoS ONE 15(3).
- 5.3 D-793, SOP, Cryogenic Grinding of Chewable Gels

6.0 Supplies

- 6.1 Chemicals
 - 6.1.1 D-(+)-Mannose and D-(-)-Ribose reference standards
 - 6.1.2 ACN
 - 6.1.3 PMP
 - 6.1.4 0.5 M HCl

Standard Operating Procedure Determination of Mannose by HPLC/UV	SOP No D-1024	Rev 1	Page 3 of 11
---	--------------------------------	------------------------	-------------------------------

6.1.5 0.5 M NaOH

6.1.6 H₃PO₄

6.1.7 K₂HPO₄

6.1.8 MeOH

6.1.9 Chloroform

6.1.10 H₂O

6.2 Glassware and Disposables

6.2.1 Volumetric glassware as required for standard and sample preparations

6.2.2 HPLC vials, 2mL with screw-cap enclosures and septa

6.2.3 Tips for adjustable pipettes

6.2.4 0.45 µm nylon or PVDF syringe filters

6.2.5 10-mL plastic syringes with luer lock fitting

6.3 Equipment

6.3.1 Suitable gradient HPLC system consisting of a pump, autosampler, column oven, and UV or DAD detector.

6.3.2 Analytical balance

6.3.3 pH meter

6.3.4 Water bath capable of maintaining temperature at 70°C

6.3.5 Wrist action shaker

Standard Operating Procedure Determination of Mannose by HPLC/UV	SOP No D-1024	Rev 1	Page 4 of 11
--	--------------------------	------------------	-------------------------

6.3.6 Adjustable pipettes

7.0 Preparation of Mobile Phase, Extraction Solvent, Standards, Internal Standards, and Samples

7.1 Mobile Phase A (15 mM K₂HPO₄ in H₂O)

7.1.1 Transfer 2.61 g of K₂HPO₄ to a suitable container.

7.1.2 Add 1000 mL of H₂O and a stir bar.

7.1.3 Begin stirring and adjust to pH 7.2 using H₃PO₄.

7.2 Mobile Phase B (ACN)

7.2.1 Use 100% ACN.

7.3 Internal Standard Solution (12.5 mg/mL D-(-)-Ribose in H₂O)

7.3.1 Accurately weigh and transfer about 625 mg of D-(-)-Ribose to a 50-mL volumetric flask.

7.3.2 Dissolve in and dilute to volume using H₂O.

7.4 Derivatization Solution (0.5 M PMP in Methanol)

7.4.1 Accurately weigh and transfer about 870 mg of PMP to a 10-mL volumetric flask.

7.4.2 Dissolve in and dilute to volume using methanol.

Standard Operating Procedure Determination of Mannose by HPLC/UV	SOP No D-1024	Rev 1	Page 5 of 11
---	--------------------------------	------------------------	-------------------------------

- 7.5 Stock Standard (2.5 mg/mL D-(+)-Mannose + 2.5 mg/mL D-(-)-Ribose)
- 7.5.1 Accurately weigh and transfer about 62.5 mg D-(+)-Mannose into a 25-mL volumetric flask.
- 7.5.2 Add 5.0 mL of Internal Standard Solution.
- 7.5.3 Dissolve in and dilute to volume using H₂O.
- 7.6 Stock Sample Preparation (~2.5 mg/mL D-(+)-Mannose + 2.5 mg/mL D-(-)-Ribose)
- 7.6.1 Specific sample testing details are provided in each products profile. If a specific testing details section is not available, then follow preparation procedure as described below, maintaining concentration within the linear range of this method.
- 7.6.2 The validated linear range of the method is 1 mg/mL – 5 mg/mL based on the concentration of mannose in the Stock Sample. The Stock Sample preparation must be within the linear range of the method.
- 7.6.3 Ensure that the sample is homogeneous prior to weighing.
- 7.6.3.1 For capsules, combine the fill material from at least 10 dosage units and homogenize in a mortar and pestle if necessary.
- 7.6.3.2 For tablets, combine at least 10 dosage units and homogenize in a mortar and pestle.
- 7.6.3.3 For chewable gels (gummies), homogenize at least 10 dosage units as outlined in D-793.
- 7.6.4 Based on the raw material manufacturer's assay value or finished product label claim, accurately weigh and transfer a sufficient amount of sample to deliver about 125 mg of D-(+)-Mannose to a 50-mL volumetric flask.

Standard Operating Procedure Determination of Mannose by HPLC/UV	SOP No D-1024	Rev 1	Page 6 of 11
---	--------------------------------	------------------------	-------------------------------

- 7.6.5 Add 5.0 mL of methanol and swirl to completely wet the sample.
- 7.6.6 Add 10.0 mL of Internal Standard Solution.
- 7.6.7 Add about 20 mL of H₂O.
- 7.6.8 Shake on a wrist-action shaker for 20 min.
- 7.6.9 Dilute to volume with H₂O, and shake vigorously.
- 7.7 Derivatization Reaction and Cleanup (Option 1)
 - 7.7.1 Transfer 1.6 mL of Stock Standard into a 15-mL glass centrifuge tube.
 - 7.7.2 Transfer 1.6 mL of Stock Sample into a separate 15-mL glass centrifuge tube.
 - 7.7.3 Add 1.6 mL of 0.5 NaOH to each of the standard/sample tubes.
 - 7.7.4 Vortex briefly to mix.
 - 7.7.5 Add 1.6 mL of Derivatization Solution to each of the standard/sample tubes.
 - 7.7.6 Vortex briefly to mix.
 - 7.7.7 Incubate in a water bath at 70°C for 90 min, shaking by hand every 30 min.
 - 7.7.8 Equilibrate to room temperature.
 - 7.7.9 Add 1.6 mL of 0.5 M HCl to each of the standard/sample tubes.
 - 7.7.10 Add 2.0 mL of chloroform to each of the standard/sample tubes.
 - 7.7.11 Vortex for about 5 sec to mix well.
 - 7.7.12 Right away vent sample to release pressure.
 - 7.7.13 Repeat steps 7.7.11 – 7.7.12 until pressure stops building up.

- 7.7.14 Centrifuge at 4000 rpm for 3 min to separate the layers.
- 7.7.15 Use a glass pipet to remove and discard the bottom (chloroform) layer.
- 7.7.16 Repeat steps 7.7.10 – 7.7.15 twice for a total of 3 washes.
- 7.8 Derivatization Reaction and Cleanup (Option 2)
 - 7.8.1 Transfer 1.6 mL of Stock Standard into a 15-mL glass centrifuge tube.
 - 7.8.2 Transfer 1.6 mL of Stock Sample into a separate 15-mL glass centrifuge tube.
 - 7.8.3 Add 1.6 mL of 0.5 NaOH to each of the standard/sample tubes.
 - 7.8.4 Vortex briefly to mix.
 - 7.8.5 Add 1.6 mL of Derivatization Solution to each of the standard/sample tubes.
 - 7.8.6 Vortex briefly to mix.
 - 7.8.7 Incubate in a water bath at 70°C for 90 min, shaking by hand every 30 min.
 - 7.8.8 Equilibrate to room temperature.
 - 7.8.9 Add 1.6 mL of 0.5 M HCl to each of the standard/sample tubes.
 - 7.8.10 Quantitatively transfer reaction mixture from standard/sample tubes to separate separatory funnels with at least 2 mL of chloroform
 - 7.8.11 Shake and vent well to release pressure, let layers separate
 - 7.8.12 Discard bottom (chloroform) layer
 - 7.8.13 Repeat steps 7.7.11 – 7.7.12 with at least 2 mL of chloroform twice more for a total of 3 washes.
 - 7.8.14 Transfer remaining (top) layer to scintillation vials

7.9 HPLC Sample Preparation

7.9.1 Transfer 1.0 mL of each reaction mixture after chloroform washing into separate 100-mL volumetric flask.

7.9.2 Dilute to volume using H₂O.

7.9.3 Filter a portion using a 0.45 µm syringe filter (nylon or PVDF), discarding the first 2-3 mL before collecting a portion in an HPLC vial for analysis.

7.10 Instrument Method Parameters

7.10.1 Column: Phenomenex Kinetex XB-C18, 4.6 x 150 mm, 2.6 µm or equivalent

7.10.2 Flow Rate: 1.0 mL/min

7.10.3 Run Time: 25 min

7.10.4 Gradient

Time (min)	%A	%B
0.0	88	12
12.0	80	20
15.0	35	65
20.0	35	65
20.1	88	12
25.0	88	12

7.10.5 Injection Volume: 5 µL

7.10.6 Column Temperature: 30 °C

7.10.7 Wavelength: 245 nm

7.10.8 Suggested 3-D Spectral Range (for Identification): 210 nm – 400 nm

Standard Operating Procedure Determination of Mannose by HPLC/UV	SOP No D-1024	Rev 1	Page 9 of 11
---	--------------------------------	------------------------	-------------------------------

7.11 Recommended Sequence

7.11.1 Make at least 1 injection of a Blank (H₂O).

7.11.2 Make 5 injections of the Working Standard.

7.11.3 Make a single injection of each Working Sample.

7.11.4 Make a single injection of the Working Standard after every six samples and at the end of the run.

7.12 System Suitability Requirements

7.12.1 The %RSD for five consecutive injections of Working Standard is NMT 3%.

7.12.2 The %RSD for all injections of Working Standard is NMT 5%.

7.12.3 No significant (>0.5%) interference is present in the Blank injection.

7.13 Column Wash and Storage

7.13.1 Wash the column at 1 mL/min with H₂O/ACN (90/10).

7.13.2 Wash the column at 1 mL/min with H₂O/ACN (50/50).

7.13.3 Store the column in H₂O/ACN (50/50).

7.14 Example Calculation

$$\% \text{ assay} = \frac{R_u}{R_s} \times \frac{Wt_{std} \times P}{V_{std}} \times \frac{SS}{Spl_{wt}} \times \frac{V_{spl}}{LA} \times 100$$

R_u Sample relative response


R_s Mean standard relative response

Wt_{std} Weight of reference standard in mg


- V_{std} Volume of the standard preparation accounting for dilutions in mL
- P Purity of the reference standard in decimal format
- SS Serving size: Weight of a single dosage unit in mg or 1 for raw materials.
- Spl_{wt} Sample weight in mg
- V_{spl} Volume of the sample preparation accounting for dilutions in mL
- LA Label amount in mg per dose or 1 for raw materials

8.0 Example Chromatograms


8.1 Blank



8.2 Standard



8.3 Sample



9.0 Revision History

Revision	Date	Description of Changes	CCR #	By
0	04/30/24	New procedure.	N/A	S. Sassman
1	10/01/25	Added use of separatory funnel as an acceptable method for derivatization step 7.8. Added clarifying steps to 7.7.	CC-25-0372	M. Autrey