	Standard Operating Procedure		SOP Number D-1025	Revision 0
	Determination of Ascorbyl Palmitate by HPLC/UV		Effective Date 05/29/24	Page Page 1 of 10
Written by/ Date SAS 05/21/24		Reviewed by/ Date CJP 05-28-24		Approved by/ Date AJS 05/29/24
Title: Analytical Development Scientist		Title: Analytical Development Scientist		Title: QC Laboratory Manager

1.0 Purpose

The purpose of this procedure is to define the method for the determination of ascorbyl palmitate in raw materials and finished products by HPLC-UV.

2.0 Scope

This procedure applies to the identification and quantification of ascorbyl palmitate in raw materials and finished products in the QC laboratory at Ion Labs.

3.0 Responsibility

- 3.1 It is the responsibility of QC Chemists who have verified their ability to execute this procedure to follow this procedure.
- 3.2 It is the responsibility of QC Laboratory Management to ensure that this procedure is being followed.
- 3.3 It is the responsibility of QC Laboratory Management and/or Analytical Development to keep this procedure aligned with current practices.

4.0 Definitions

- 4.1 QC – Quality control
- 4.2 HPLC-UV – High Performance Liquid Chromatography with Ultraviolet Detection
- 4.3 ACN – Acetonitrile

Standard Operating Procedure Determination of Ascorbyl Palmitate by HPLC/UV	SOP No D-1025	Rev 0	Page 2 of 10
--	--------------------------------	------------------------	-------------------------------

4.4 **H₃PO₄** – Phosphoric acid

4.5 **H₂O** – Deionized water (>18MΩ·cm)

5.0 References

5.1 PRTCL-24-0038, Protocol, Validation of a Method for the Determination of Ascorbyl Palmitate by HPLC/UV

5.2 D-903, SOP, Conversion Factors Used in Analytical Determinations and New Product Formulations

6.0 Supplies

6.1 Chemicals

6.1.1 Ascorbyl palmitate reference standard

6.1.2 ACN

6.1.3 H₃PO₄

6.1.4 Methanol

6.1.5 Citric acid

6.1.6 Ascorbic acid

6.2 Glassware and Disposables

6.2.1 Volumetric glassware as required for standard and sample preparations

6.2.2 HPLC vials, 2mL with screw-cap enclosures and septa

6.2.3 Tips for adjustable pipettes

6.2.4 0.45 μm nylon or PVDF syringe filters

Standard Operating Procedure Determination of Ascorbyl Palmitate by HPLC/UV	SOP No D-1025	Rev 0	Page 3 of 10
---	--------------------------	------------------	-------------------------

6.2.5 10-mL plastic syringes with luer lock fitting

6.3 Equipment

6.3.1 Suitable gradient HPLC system consisting of a pump, autosampler, column compartment and UV-Vis detector with a chromatographic data handling system

6.3.2 Analytical balance

6.3.3 Micro balance

6.3.4 Wrist action shaker

6.3.5 Adjustable pipettes

7.0 Preparation of Mobile Phase, Diluent, Standards, and Samples

7.1 Mobile Phase A (0.1% H₃PO₄)

7.1.1 Transfer 1000 mL of H₂O to a suitable container.

7.1.2 Add 1.0 mL of H₃PO₄, and mix well.

7.1.3 Scale as necessary. Mobile phase A has an expiry of one month.

7.2 Mobile Phase B

7.2.1 Use 100% ACN.

7.2.2 Mobile phase B has an expiry of three months.

7.3 Diluent A (1 g/L citric acid + 1g/L ascorbic acid in methanol)

7.3.1 Transfer 1.0 g of citric acid to a suitable container.

7.3.2 Add 1.0 g of ascorbic acid.

7.3.3 Add 1000 mL of methanol.

<p style="text-align: center;">Standard Operating Procedure Determination of Ascorbyl Palmitate by HPLC/UV</p>	<p style="text-align: center;">SOP No D-1025</p>	<p style="text-align: center;">Rev 0</p>	<p style="text-align: center;">Page 4 of 10</p>
--	---	---	--

- 7.3.4 Sonicate or stir until completely dissolved.
- 7.3.5 Equilibrate to room temperature before use.
- 7.3.6 Store at 4°C. The diluent has an expiry of one month.
- 7.4 Diluent B (1 g/L citric acid + 1g/L ascorbic acid in DMSO/methanol 50/50)
 - 7.4.1 Diluent B is used for chewable gels (gummies) only.
 - 7.4.2 Transfer 1.0 g of citric acid to a suitable container.
 - 7.4.3 Add 1.0 g of ascorbic acid.
 - 7.4.4 Add 500 mL of methanol.
 - 7.4.5 Add 500 mL of DMSO
 - 7.4.6 Sonicate or stir until completely dissolved.
 - 7.4.7 Equilibrate to room temperature before use.
 - 7.4.8 Store at 4°C. The diluent has an expiry of one month.
- 7.5 Stock Standard (~500 mcg/mL)
 - 7.5.1 Accurately weigh and transfer about 25 mg of reference standard into a 50-mL low-actinic (red) volumetric flask.
 - 7.5.2 Dissolve in and dilute to volume using Diluent A.
- 7.6 Working Standard (~100 mcg/mL)
 - 7.6.1 Transfer 10.0 mL of Stock Standard to a 50-mL low-actinic (red) volumetric flask.
 - 7.6.2 Dilute to volume with Diluent A.

<p style="text-align: center;">Standard Operating Procedure Determination of Ascorbyl Palmitate by HPLC/UV</p>	<p style="text-align: center;">SOP No D-1025</p>	<p style="text-align: center;">Rev 0</p>	<p style="text-align: center;">Page 5 of 10</p>
--	---	---	--

7.7 Sample Preparation

- 7.7.1 Specific sample testing details are provided in each products profile. If a specific testing details section is not available, then follow preparation procedure as described below, maintaining concentration within the linear range of this method.
- 7.7.2 The linear range is 5 – 500 mcg/mL of ascorbyl palmitate.
- 7.7.3 Ensure that the sample is homogeneous prior to weighing.
 - 7.7.3.1 For capsules, combine the fill material from at least 10 dosage units and homogenize in a mortar and pestle if necessary.
 - 7.7.3.2 For tablets, combine at least 10 dosage units and homogenize in a mortar and pestle.
 - 7.7.3.3 For chewable gels (gummies), homogenize at least 10 dosage units as outlined in D-793.
- 7.7.4 For raw materials: Weigh no less than 25 mg into a suitably sized volumetric flask of no less than 25 mL volume to generate an analyte concentration that is within the validated linearity range. Fill the flask to about 65% of the calculated volume with Diluent A, shake mechanically for 20 minutes, and dilute to volume using Diluent A.
- 7.7.5 For solid and liquid dose finished products: Based on the label claim and fill weight (capsules), serving size (powders and liquids) or tablet weight per dose, weigh no less than 100 mg of the pooled dosages into a suitably sized volumetric flask of no less than 25 mL to generate an analyte concentration that is within the validated linear range. Fill the flask to about 65% of the calculated volume with Diluent A, shake mechanically for 20 minutes, and dilute to volume using Diluent A.

Standard Operating Procedure Determination of Ascorbyl Palmitate by HPLC/UV	SOP No D-1025	Rev 0	Page 6 of 10
--	--------------------------------	------------------------	-------------------------------

- 7.7.6 For chewable gels (gummies): Homogenize at least 10 dosage units according to the procedure outlined in D-793 Cryogenic Grinding of Chewable Gels. Quickly weigh no less than 200 mg of the homogenized sample into a suitably sized beaker. Add a volume of Diluent B equivalent to 50% of the desired flask volume, add a stir bar, protect from light, and stir for 30 min. Transfer the resulting solution to a volumetric flask of size suitable to generate an analyte concentration that is within the validated linear range. Use several small portions of Diluent B to rinse any remaining residue from the beaker into the volumetric flask ensuring complete transfer, then sonicate for 10 min. Equilibrate to room temperature, and dilute to volume using Diluent B.
- 7.7.7 To manage large volumes, the sample can be initially dissolved in a smaller volume and a portion further diluted using Diluent A to bring the analyte concentration into the linear range. Dilutions can be made using volumetric glassware and/or adjustable pipettes. Dilutions can be prepared in HPLC vials.
- 7.7.8 Filter through a 0.45 μm membrane discarding the first 2 – 3 mL before collecting a portion for analysis. Alternatively, centrifuge an aliquot of the final sample at 10,000 rpm for 5 min to remove particulates.

7.8 Instrument Method Parameters

- 7.8.1 Column: Ascentis Express C8, 2.7 μm , 4.6 mm x 100 mm
- 7.8.2 Mobile phase: Isocratic A/B (33/67)
- 7.8.3 Flow Rate: 1.0 mL/min
- 7.8.4 Injection volume: 10 μL
- 7.8.5 Column temperature: 30 $^{\circ}\text{C}$
- 7.8.6 Wavelength: 242 nm
- 7.8.7 Recommended Spectral Range (for Identification): 200 – 400 nm

Standard Operating Procedure Determination of Ascorbyl Palmitate by HPLC/UV	SOP No D-1025	Rev 0	Page 7 of 10
--	--------------------------------	------------------------	-------------------------------

7.9 Recommended Sequence

7.9.1 Make at least 2 injections of Diluent.

7.9.2 Make 5 injections of the Working Standard.

7.9.3 Make a single injection of each Working Sample.

7.9.4 Make a single injection of the Working Standard after every six samples and at the end of the run.

7.10 System Suitability Requirements

7.10.1 The %RSD for five consecutive injections of Working Standard is NMT 2%.

7.10.2 The %RSD for all injections of Working Standard is NMT 3%.

7.10.3 No significant (>0.5%) interferences are present in the diluent injection.

7.11 Column Cleaning

7.11.1 It is recommended to clean the column after every run.

7.11.2 Clean the column at 0.2 mL/min with Mobile Phase A / Mobile Phase B (20/80) for at least 30 min.

7.12 Column Storage

7.12.1 Wash and store the column in ACN/H₂O (50/50).

8.0 Example Calculation

$$\% \text{ assay} = \frac{R_u}{R_s} \times \frac{W_{t_{\text{std}}} \times P}{V_{\text{std}}} \times \frac{SS}{Spl_{wt}} \times \frac{V_{\text{spl}}}{LA} \times \text{MULT} \times 100$$


R_u Sample peak area

R_s Mean standard peak area


Wt _{std}	Weight of reference standard in mg
V _{std}	Volume of the standard preparation accounting for dilutions in mL
P	Purity of the reference standard in decimal format
SS	Serving size: Weight of a single dosage unit in mg or 1 for raw materials.
Spl _{wt}	Sample weight in mg
V _{spl}	Volume of the sample preparation accounting for dilutions in mL
LA	Label amount in mg per dose or 1 for raw materials
MULT	Multiplier (from SOP D-903 Conversion Factors Used in Analytical Determinations and New Product Formulations)

9.0 Example Chromatograms


9.1 Blank



9.2 Standard



9.3 Sample



Standard Operating Procedure Determination of Ascorbyl Palmitate by HPLC/UV	SOP No D-1025	Rev 0	Page 10 of 10
--	--------------------------------	------------------------	--------------------------------

10.0 Revision History

Revision	Date	Description of Changes	CCR #	By
0	05/02/24	New procedure.	N/A	S. Sassman