

	<b>Standard Operating Procedure</b> <b>Ascorbic Acid Determination by HPLC</b> <b>with UV/Vis Spectroscopy</b>	<b>SOP Number</b> <b>D-701</b>	<b>Revision</b> <b>1</b>
		<b>Effective Date</b> 01/03/23	<b>Page</b> <b>Page 1 of 9</b>
<b>Written by/ Date</b>  12/20/22	<b>Reviewed by/ Date</b> SAS 12/20/22	<b>Approved by/ Date</b>  12/20/22	
<b>Title: Quality Control</b> <b>Director</b>	<b>Title: Analytical Development</b> <b>Scientist</b>	<b>Title: Quality Assurance</b> <b>Director</b>	

## 1.0 Purpose

This document describes the analytical procedure for the determination of Ascorbic Acid (AA) in raw materials and finished products.

## 2.0 Scope

This procedure applies to the identification and quantification of AA in raw materials and finished products. This method was validated under Protocol PRTCL-20-0002.

## 3.0 Responsibility

- 3.1 It is the responsibility of QC and Analytical chemists who have verified their ability to execute this procedure to follow this procedure.
- 3.2 It is the responsibility of QC Laboratory Management to implement this procedure and to ensure that the procedure is being followed.
- 3.3 It is the responsibility of QC Laboratory Management and AD Personnel to keep this procedure current with the associated monographs and laboratory practices.

## 4.0 Definitions

- 4.1 **QC** – Quality Control
- 4.2 **AD** – Analytical Development
- 4.3 **AA** – Ascorbic Acid
- 4.4 **ACN** – Acetonitrile

- 4.5 **TFA** – Trifluoroacetic Acid
- 4.6 **EDTA** – Ethylenediaminetetraacetic Acid
- 4.7 **ACS** – American Chemical Society
- 4.8 **HPLC** – High Performance Liquid Chromatography
- 4.9 **UV/Vis** – Ultraviolet & Visible Electromagnetic Spectra

## **5.0 References**

- 5.1 PRTCL-20-0002, Protocol, Ascorbic Acid Determination by HPLC using UV/Vis Spectroscopy
- 5.2 D-793, SOP, Cryogenic Grinding of Chewable Gels

## **6.0 Supplies**

- 6.1 Chemicals – All reagents are ACS grade or better.
  - 6.1.1 Milli-Q Water
  - 6.1.2 ACN
  - 6.1.3 TFA
  - 6.1.4 Ascorbic Acid Reference Standard
  - 6.1.5 EDTA Disodium Dihydrate
  - 6.1.6 Sodium Phosphate Monobasic Dihydrate
  - 6.1.7 Phosphoric Acid
- 6.2 Supplies and Glassware
  - 6.2.1 HPLC vials, 12mm X 32mm with screw cap enclosures w/ septa

6.2.2 Volumetric glassware and/or adjustable pipettes and tips

6.2.3 Weigh paper or funnels

6.2.4 10ml Syringes with 17mm x 0.45u Nylon or Glass Fiber Syringe Filters (Note: It is recommended that chewable gel samples be filtered using glass fiber syringe filters.)

### 6.3 Equipment

6.3.1 Suitable gradient HPLC system consisting of a pump, autosampler, column oven and UV detector with a chromatographic data handling system

6.3.2 Analytical Balance

6.3.3 Wrist Action Shaker

## 7.0 Procedure

### 7.1 Mobile Phase & Diluent Preparation

#### 7.1.1 Mobile Phase

7.1.1.1 Mobile Phase A: Add 200  $\mu$ L of TFA to 1000 mL of water and mix well.

7.1.1.2 Mobile Phase B: ACN

#### 7.1.2 Extraction Solvent / Diluent

7.1.2.1 Dissolve 2.34g Sodium Phosphate Monobasic Dihydrate and 0.56g ETDA Disodium Dihydrate in 1000ml water then pH to 3.0 with Phosphoric Acid.

7.1.3 Preparations may be scaled as necessary

### 7.2 Standard Prep

- 7.2.1 Accurately weigh and transfer about 45 mg of AA reference standard into a 100-mL volumetric flask. Add 50mL of Diluent and briefly sonicate until all AA is dissolved. (**Caution:** Perform this step quickly as AA is air sensitive.)
- 7.2.2 Dilute to volume with Diluent and mix well – this is the AA Stock. Dilute the AA Stock 1:10 with Diluent – this is the AA Working Standard.
- 7.2.3 Alternative standard preparations are acceptable as long as the preparations are within the linear range of this method, 0.01856 – 0.09280 mg/mL

### 7.3 Sample Preparation

- 7.3.1 Specific sample testing details are provided in each products profile. If a specific testing details section is not available, then follow preparation procedure as described below, maintaining concentration within the linear range of this method.
- 7.3.2 The validated range for the analytical method is 0.01856 – 0.09280 mg/mL.
- 7.3.3 For finished products, extract sufficient sample with Diluent in order to generate an AA concentration that is within the validated linearity range. (**Caution:** Some finished products contain alkaline excipients in quantities capable of overwhelming the buffer capacity of the Diluent. Make sure that extraction sample solutions never exceed a pH of 3.0.) When analyzing gummies, powder and dissolve dosage forms as detailed in D-793 Cryogenic Grinding of Chewable Gels.
- 7.3.4 Prepare raw materials like standards. (However, be sure to consult the specification for expected potency, as raw material samples may not be 100%.)
- 7.3.5 Samples can be dissolved in Diluent at any volume starting from 100mL. The volume chosen must be in the solubility range of AA (validated at approximately 0.4 mg/ml). To manage large volumes, the sample can be initially dissolved in a smaller volume that is within the solubility range and a portion further diluted to

bring the AA concentration into the linear range.

7.3.6 Fill the flask to about 50% of the calculated volume with Diluent and shake mechanically for 10 minutes. (Note: finished products that are liquids, along with chewable gel samples that have been previously powdered and dissolved in Diluent as per D-793 Cryogenic Grinding of Chewable Gels do not require mechanical shaking.) QS to volume with Diluent.

7.3.7 Perform further dilutions as required using Diluent. Filter a 5ml aliquot for analysis, discarding the first 3-4ml of filtrate.

#### 7.4 HPLC Parameters

7.4.1 Column: Phenomenex Kinetex XB-C18, 4.6 x 150mm, 2.6µm (Or Equivalent)

7.4.2 Column Temperature: 30°C

7.4.3 Flow rate: 0.7 mL/min

7.4.4 Mobile Phase Gradient:

Time, min	% A	% B
0.0	100	0
1.0	100	0
6.0	95	5
6.1	100	0
12.0	100	0

7.4.5 Wavelength: 254 nm

7.4.6 Injection Volume: 20 µL

7.4.7 Run Time: 12 minutes (Note: Some finished products may require intermittent blank or flush injections in order to remove more highly retentive formulation constituents that could interfere with subsequent injections.)

7.4.8 Recommended 3-D Spectral Range (for Identification) - 210nm to 350nm

7.5 Recommended Sequence

7.5.1 Make at least 2 injections of the Diluent.

7.5.2 Make at least five (5) injections of AA Working Standard.

7.5.3 Make a single injection of each Sample Preparation.

7.5.4 Make a single injection of the Standard Solution after every ten (10) sample injections and/or at the end of a run.

7.6 System Suitability Requirements

7.6.1 The %RSD of five (5) consecutive standard injections is NMT 2.0%

7.6.2 The %RSD of all standard injections is NMT 3.0%.

7.6.3 If present, any interference in the diluent should be subtracted out of the sample and standard peak areas.

7.7 Example calculations for determining finished product % label claim or raw material % assay:

$$7.7.1 \quad \% \text{ AA} = \frac{R_u}{R_s} \times \frac{Wt_{std} \times P}{V_{std}} \times \frac{SS}{SA} \times \frac{V_{spl}}{LA} \times 100$$

$R_u$             Sample peak area

$R_s$             Mean (n=5) standard peak area

$Wt_{std}$         Weight of the reference standard (mg)

$V_{std}$         Volume of the standard preparation accounting for dilutions in mL

$P$              Purity of the reference standard in decimal format

$SA$             Sample amount (mg)

$V_{spl}$  Volume of the sample preparation accounting for dilutions in mL

$SS$  Serving size: Average weight of ten dosage units for tablets and chewable gels, fill weight for capsules, mass of a single serving for powders, volume of a single serving from the theoretical formula in mL for liquids, or 1 for raw materials.


$LA$  Label amount of AA (use 1 for raw materials)

## 7.8 System Wash, Column Wash and Column Storage


7.8.1 Wash and store the column in 50:50 ACN / Water.

## 8.0 Chromatograms


### 8.1 Typical Diluent Chromatogram




### 8.2 Typical Working Standard Chromatogram



### 8.3 Typical Raw Material Chromatogram



### 8.4 Typical Finished Product Chromatogram



### 9.0 Revision History

Revision	Date	Description of Changes	CCR #	By
0	02/12/20	New	N/A	C. Perry
1	12/20/22	Add sample test details information. Minor edits.	CC-22-0475	J. Sassman