	Standard Operating Procedure Corosolic Acid Determination by HPLC using UV/VIS Spectroscopy		SOP Number D-766	Revision 1
			Effective Date 05/13/22	Page Page 1 of 7
Written by/ Date SAS 05/02/22		Reviewed by/ Date Jm 05/02/22		Approved by/ Date SS 05/02/22
Title: Analytical Development Scientist		Title: Analytical Development Manager		Title: QC Laboratory Director

1.0 Purpose

The purpose of this procedure is to define the method for the quantitation and/or identification of corosolic acid in raw materials and finished product dietary supplements using HPLC and UV/VIS spectrophotometry.

2.0 Scope

This procedure applies to the quantification and identification of corosolic acid in raw materials and finished products. Corosolic acid is a decent chromophore and was measured at 210, other wavelengths can be used to maximize signal to noise.

3.0 Responsibility

- 3.1 It is the responsibility of QC Chemists to follow this procedure.
- 3.2 It is the responsibility of QC Laboratory Management to ensure that this procedure is being followed.
- 3.3 It is the responsibility of QC Laboratory and Analytical Development Management to keep this procedure aligned with current practices.

4.0 Definitions

- 4.1 **UV/VIS** – Ultraviolet and Visible Electromagnetic Spectrums
- 4.2 **H₃PO₄** – Phosphoric acid
- 4.3 **ACN** – Acetonitrile
- 4.4 **CofA** – Certificate of Analysis
- 4.5 **H₂O** – Water ($\geq 18.2 \text{ M}\Omega\cdot\text{cm}$)

Standard Operating Procedure Corosolic Acid Determination by HPLC using UV/VIS Spectroscopy	SOP No D-766	Rev 1	Page 2 of 7
--	-------------------------------	------------------------	------------------------------

- 4.6 **Corosolic acid** – (1*S*,2*R*,4*aS*,6*aR*,6*aS*,6*bR*,8*aR*,10*R*,11*R*,12*aR*,14*bS*)-10,11-Dihydroxy-1,2,6*a*,6*b*,9,9,12*a*-heptamethyl-2,3,4,5,6,6*a*,7,8,8*a*,10,11,12,13,14*b*-tetradecahydro-1*H*-picene-4*a*-carboxylic acid

5.0 References

- 5.1 MV-LAB-18-206, Protocol, Corosolic Acid Determination by HPLC using UV/Vis Spectroscopy
- 5.2 Vijaykumar, Katta, et al. Int. J. Appl. Sci. Eng. 2006. Vol 4. Pg. 103-114.
- 5.3 USP 43 Monograph for Banaba Leaf

6.0 Supplies

- 6.1 Chemicals: All reagents are HPLC grade or better.
- 6.1.1 H₂O
- 6.1.2 ACN
- 6.1.3 H₃PO₄
- 6.1.4 Corosolic acid reference standard
- 6.1.5 Methanol
- 6.2 Glassware
- 6.2.1 HPLC vials, 12mm x 32mm with screw cap enclosures with septa
- 6.2.2 Scintillation Vials
- 6.2.3 1L Mobile Phase Container
- 6.2.4 50mL Volumetric Flask
- 6.2.5 100mL Volumetric Flask
- 6.3 Disposables
- 6.3.1 10mL Pipette Tips
- 6.3.2 1mL Pipette Tips
- 6.3.3 200µL Pipette Tips

Standard Operating Procedure Corosolic Acid Determination by HPLC using UV/VIS Spectroscopy	SOP No D-766	Rev 1	Page 3 of 7
--	-------------------------------	------------------------	------------------------------

6.3.4 1.5mL microfuge tubes

6.3.5 16mL Test Tubes

6.3.6 Disposable Plastic Luer Lock Syringe – 3mL, 6mL, or 10mL

6.3.7 Nylon Syringe Filters, 0.45µm

6.3.8 Weigh paper

6.4 Equipment

6.4.1 Suitable gradient HPLC system consisting of a pump, autosampler, column oven and UV detector with a chromatographic data handling system

6.4.2 Analytical Balance

6.4.3 Ultrasonic bath

6.4.4 Vortex

6.4.5 Stir Plate

6.4.6 Eppendorf Centrifuge

6.4.7 10mL Pipette

6.4.8 1mL Pipette

6.4.9 200µL Pipette

7.0 Preparation of Mobile Phase, Dissolution Buffer, Samples, and Standards

7.1 Mobile Phase A – 0.1% Phosphoric acid in H₂O

Prepared by adding 1.0mL of H₃PO₄ to 1000mL of H₂O (scale as necessary).

7.2 Mobile Phase B – 0.1% Phosphoric acid in Acetonitrile

Prepared by adding 1.0mL of H₃PO₄ to 1000mL of ACN (scale as necessary).

7.3 Diluent– Methanol

7.4 Standard Preparation

7.4.1 The linear range of the method is 0.01 mg/mL – 0.1 mg/mL. All final standard and sample preparations must be within the linear range.

Standard Operating Procedure Corosolic Acid Determination by HPLC using UV/VIS Spectroscopy	SOP No D-766	Rev 1	Page 4 of 7
--	-------------------------------	------------------------	------------------------------

7.4.2 Use the actual purity from the CofA or the standard certification for the reference material for calculations.

7.4.3 The standard is prepared by weighing no less than the minimum weight of the analytical balance, then bringing up to the final volume in an appropriate volumetric flask using Diluent, then sonicating for 10 minutes.

7.4.4 Dilutions are prepared using Diluent. Dilutions can be made using volumetric flasks or using 10mL, 1mL, and 200µL variable pipettes. Specific standard concentrations will approximate the concentration expected to be found in the product being tested based on the sample dilution and calculated from the label. Dilutions can be prepared in HPLC vials.

7.5 Sample Preparation

7.5.1 Samples can be dissolved in Diluent at any volume starting from 50mL and any weight greater than the minimum weight of the analytical balance.

7.5.2 The sample is suspended in the final volume and put in the sonicator bath for at least 10 minutes with occasional shaking to ensure large chunks are dispersed.

7.5.3 Before injection, insoluble matter should be removed via filtration using a nylon syringe filter. Discard at least the first 0.5mL of filtrate before collecting a sample for analysis. Dilute filtrate as needed then add 1mL of the final dilution to an HPLC vial for analysis.

7.5.3.1 Alternatively, samples and standards can also be centrifuged at 6,000 RPM for at least 200 seconds in an Eppendorf centrifuge to pellet insoluble matter.

7.5.4 For raw materials or finished products being analyzed for the first time using this method, in-process verification is required to demonstrate spectral purity and extraction efficiency before the method can be implemented.

8.0 Test Conditions

8.1 Gradient (Isocratic)

Time	%A	%B
------	----	----

Standard Operating Procedure Corosolic Acid Determination by HPLC using UV/VIS Spectroscopy	SOP No D-766	Rev 1	Page 5 of 7
--	-------------------------------	------------------------	------------------------------

0.00 25 75

20.0 25 75

8.2 Column – Phenomenex Luna, C18, 5um, 100Å, LC column, 250mm x 4.6mm, or equivalent.

8.3 Flow Rate – 1.0mL/min

8.4 UV Detection – 210 nm

8.5 3D Spectral Range – 200 nm – 350 nm

8.6 Injection Volume - 20µL

8.7 Retention Time – about 11 minutes

8.8 Column Temperature – Not controlled

8.9 Recommended Sequence

8.9.1 Make at least two injections of the diluent.

8.9.2 Make at least two injections of the diluent.

8.9.3 Make five injections of Standard Solution.

8.9.4 Make a single injection of each Sample Preparation.

8.9.5 Make a single injection of the Standard Solution after every six samples and at the end of the run.

8.10 System Suitability Requirements

8.10.1 The %RSD of the first five standard injections is NMT 5.0%.

8.10.2 The %RSD of all injections is NMT 5%.

8.10.3 Spectral match over the range 200 nm – 350 nm is NLT 900.

8.10.4 The retention time of the sample is within 0.3 min of the standard.

8.11 Column Wash and Storage

8.11.1 Wash the column with H₂O:ACN (25:75) at 1 mL/min for at least 15 min.

8.11.2 Store the column with H₂O:ACN (25:75).

Standard Operating Procedure Corosolic Acid Determination by HPLC using UV/VIS Spectroscopy	SOP No D-766	Rev 1	Page 6 of 7
--	-------------------------------	------------------------	------------------------------

9.0 Calculations

$$9.1 \quad \% \text{ assay} = \frac{R_u}{R_s} \times \frac{W_{t_{std}} \times P}{V_{std}} \times \frac{V_{spl}}{SA} \times \frac{SS}{LA} \times 100$$

R_u Sample peak area

R_s Mean standard peak area

$W_{t_{std}}$ Weight of reference standard in mg

V_{std} Volume of the standard preparation accounting for dilutions in mL

P Purity of the reference standard in decimal format

SA Sample amount in mg (solids) or mL (liquids)


V_{spl} Volume of the sample preparation accounting for dilutions in mL

SS Serving size: Weight of a single dosage unit in mg for tablets and capsules, volume of a single serving from the theoretical formula in mL for liquids, or 1 for raw materials.


LA Label amount in mg per dose or 1 for raw materials

10.0 Example Chromatography


10.1 Blank



10.2 Standard



10.3 Sample



11.0 Revision History

Revision	Date	Description of Changes	CCR #	By
0	01/02/19	New	N/A	J. Maignan
1	04/15/22	Update to reflect current practices and for clarity, add reference to method validation, add linear range, add recommended sequence, add system suitability requirements, add example chromatography.	CC-22-0179	S. Sassman